


Anti-Efnb1 Rabbit Polyclonal Antibody For Wb, Elisa - P98172

Item Number: CM0014096



Introduction

High-quality rabbit polyclonal antibody targeting human EFN1 (Ephrin-B1), suitable for Western blot and ELISA; cross-reacts with mouse. Ideal for investigating neuronal, vascular, and epithelial development.

[Learn More](#)

Parameter	Value
Product Name	EFNB1 Rabbit pAb
Remarks/Alias	CFND; CFNS; EFB1; EFL3; EPLG2; Elk-L; LERK2; EFN1
Species	Human
Gene ID	1947
Immunogen	Recombinant protein/Recombinant fusion protein containing a sequence corresponding to amino acids 259-346 of human EFN1 (NP_004420.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P98172
Protein Weight	38kDa
Shipping	Ice bag