

Anti-Efs Rabbit Polyclonal Antibody For Wb And Elisa - O43281

Item Number: CM0024666

Introduction

Rabbit polyclonal antibody targeting EFS (Embryonal Fyn-associated substrate, CASS3) validated for Western blot and ELISA, with cross-reactivity to mouse. Ideal for cell adhesion signaling research.

[Learn More](#)

Parameter	Value
Product Name	EFS Rabbit pAb
Remarks/Alias	SIN; CAS3; EFS1; EFS2; HEFS; CASS3; EFS
Species	Human
GeneID (Human)	10278
GeneID	10278
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-190 of human EFS (NP_115835.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O43281
Protein Weight	59kDa
Shipping	Ice bag