

Anti-Eftud2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q15029

Item Number: CM0029626



Introduction

Rabbit polyclonal antibody targeting human EFTUD2 (U5-116 kDa, SNU114), validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with mouse and rat, suitable for spliceosome and mRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	EFTUD2 Rabbit pAb
Remarks/Alias	MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-116KD; EFTUD2
Species	Human
Gene ID (Human)	9343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human EFTUD2 (NP_004238.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15029
Protein Weight	109kDa
Shipping	Ice bag