

Anti-Eftud2 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q15029

Item Number: CM0004978



Introduction

Validated anti-EFTUD2 rabbit monoclonal antibody for detection of human EFTUD2 (spliceosome component, nucleus) in WB, IHC-P, IF/ICC, IP, ELISA; cross-reacts with mouse, rat. Ideal for pre-mRNA splicing studies.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EFTUD2 Rabbit mAb
Remarks/Alias	MFDm; MFDGA; Snr114; Snrp116; SNRNP116; U5-116KD
Species	Human
GeneID (Human)	9343
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-205 of human EFTUD2 (NP_004238.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15029
Protein Weight	109kDa
Shipping	Ice bag