

Anti-Eg5/Kif11 Rabbit Monoclonal Antibody For Wb & Elisa - P52732

Item Number: CM0010643



Introduction

Anti-EG5/KIF11 rabbit monoclonal antibody suitable for WB and ELISA applications. Targets human Eg5 kinesin involving in mitosis and intracellular transport. 119 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	EG5/KIF11 Rabbit mAb
Remarks/Alias	EG5; HKSP; KNSL1; MCLMR; TRIP5; EG5/KIF11
Species	Human
GeneID (Human)	3832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 900-1056 of human EG5/KIF11 (P52732)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P52732
Protein Weight	119kDa
Shipping	Ice bag