

## Anti-Eg5/Kif11 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52732

Item Number: CM0027467



### Introduction

CamelBio's EG5/KIF11 Rabbit pAb validates in WB, IF/ICC, ELISA, targeting human EG5 (119kDa, P52732). A kinesin motor protein critical for spindle bipolarity, mitosis, and Golgi secretion.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | EG5/KIF11 Rabbit pAb                                                                                                        |
| Remarks/Alias    | EG5; HKSP; KNSL1; MCLMR; TRIP5; EG5/KIF11                                                                                   |
| Species          | Human                                                                                                                       |
| Gene ID (Human)  | 3832                                                                                                                        |
| Gene ID          | 3832                                                                                                                        |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 797-1056 of human EG5/EG5/KIF11 (NP_004514.2) |
| Source           | Rabbit                                                                                                                      |
| Category         | Polyclonal Antibodies                                                                                                       |
| Application      | WB, IF/ICC, ELISA                                                                                                           |
| Cross Reactivity | Human                                                                                                                       |
| SWISS            | P52732                                                                                                                      |
| Protein Weight   | 119kDa                                                                                                                      |
| Shipping         | Ice bag                                                                                                                     |