

Anti-Egf Monoclonal Antibody For Wb - P01132

Item Number: CM0012518



Introduction

High-quality rabbit monoclonal antibody against mouse EGF, a growth factor that stimulates cell proliferation and magnesium reabsorption. Validated for Western blot and ELISA. Perfect for research in cell signaling, renal function, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	EGF Rabbit mAb
Remarks/Alias	EGF
Species	Mouse
Gene ID (Human)	1950
Gene ID	13645
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 977-1029 of mouse EGF (NP_034243.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P01132
Protein Weight	133kDa
Shipping	Ice bag