

Anti-Egf Polyclonal Antibody For Wb, If-P, Elisa - P01133

Item Number: CM0014584



Introduction

Rabbit polyclonal antibody against human EGF (Pro-epidermal growth factor) for WB, IF-P, and ELISA applications. Cross-reacts with mouse and rat EGF. Suitable for growth factor and tissue regeneration research.

[Learn More](#)

Parameter	Value
Product Name	EGF Rabbit pAb
Remarks/Alias	URG; HOMG4; EGF
Species	Human
Gene ID (Human)	1950
Gene ID	1950
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 851-950 of human EGF (NP_001954.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01133
Protein Weight	134kDa
Shipping	Ice bag