

Anti-Egf Rabbit Polyclonal Antibody For Wb, Db, Elisa - P01133

Item Number: CM0030425



Introduction

Rabbit polyclonal antibody against human EGF (Pro-epidermal growth factor). Validated for WB, DB, ELISA; cross-reacts with mouse.

Recombinant immunogen. For research use in growth factor signaling and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	EGF Rabbit pAb
Remarks/Alias	URG; HOMG4; EGF
Species	Human
GeneID (human)	1950
GeneID	1950
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 971-1023 of human EGF (NP_001954.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01133
Protein Weight	134kDa
Shipping	Ice bag