

Anti-Egfl7 Rabbit Polyclonal Antibody For Wb - Q9Uhf1

Item Number: CM0029648

Introduction

Anti-EGFL7 rabbit polyclonal antibody validated for Western blot (WB) and ELISA, cross-reacting with mouse and rat. Suitable for angiogenesis and vascular biology research. Target protein weight ~30 kDa.

[Learn More](#)

Parameter	Value
Product Name	EGFL7 Rabbit pAb
Remarks/Alias	NEU1; ZNEU1; VE-STATIN; EGFL7
Species	Human
GeneID (Human)	51162
GeneID	51162
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-100 of human EGFL7 (NP_958854.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UHF1
Protein Weight	30kDa
Shipping	Ice bag