


Anti-Egfr (L858R) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P00533

Item Number: CM0010750



Introduction

EGFR (L858R) rabbit monoclonal antibody specific for the L858R mutant of human EGFR. Suitable for Western blot, IHC-paraffin, and ELISA; cross-reacts with mouse. Ideal for studying oncogenic signaling and tyrosine kinase inhibitor resistance.

[Learn More](#)

Parameter	Value
Product Name	EGFR (L858R) Rabbit mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; EGFR (L858R)
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 800-900 of human EGFR (L858R) (P00533).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag