

Anti-Egfr Monoclonal Antibody For Wb, If-P, Elisa - P00533

Item Number: CM0004516



Introduction

KD-validated anti-EGFR rabbit monoclonal antibody for Western blot, IF-paraffin, and ELISA. Specific to human, mouse, and rat EGFR. Recombinant immunogen-based, ideal for cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EGFR Rabbit mAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; EGFR
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag