

Anti-Egfr (Phospho Y1068) Rabbit Monoclonal Antibody For Wb And Elisa - P00533

Item Number: CM0000563



Introduction

Phospho-EGFR-Y1068 Rabbit mAb is a rabbit monoclonal antibody specific for EGFR phosphorylated at tyrosine 1068. Validated for WB and ELISA, it cross-reacts with human and mouse, making it ideal for EGFR signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-EGFR-Y1068 Rabbit mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; Phospho-EGFR-Y1068
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y1068 of human EGFR (P00533).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag