

Anti-Egfr Rabbit Monoclonal Antibody For Fc - P00533

Item Number: CM0011069



Introduction

Rabbit anti-Human EGFR mAb detects human EGFR (P00533) by flow cytometry. This rabbit monoclonal antibody targets the extracellular domain of EGFR, a 134 kDa receptor tyrosine kinase pivotal in cell signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human EGFR mAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag