


Anti-Egfr Rabbit Monoclonal Antibody For Wb, Elisa - P00533

Item Number: CM0004535



Introduction

KD validated anti-EGFR rabbit monoclonal antibody for Western blot and ELISA. Specific for human EGFR (134 kDa), suitable for cancer signaling and receptor research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EGFR Rabbit mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2; FR
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 363-463 of human EGFR. (NP_005219.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag