

Anti-Egfr Rabbit Polyclonal Antibody [Kd Validated] For Wb, Ihc-P, If/Icc, Elisa - P00533

Item Number: CM0016229



Introduction

High-quality rabbit polyclonal antibody targeting human EGFR (UniProt P00533). Validated for WB, IHC-P, IF/ICC, and ELISA. Ideal for investigating EGFR signaling in cancer, cell migration, and viral entry research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EGFR Rabbit pAb
Remarks/Alias	ERBB; ERBP; HER1; mENA; ERBB1; PIG61; NISBD2; FR
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1021-1210 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P00533
Protein Weight	134kDa
Shipping	Ice bag