

Anti-Egln1/Egln2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Gzt9/Q96Ks0

Item Number: CM0013641



Introduction

Rabbit polyclonal antibody targeting human EGLN1/EGLN2 (PHD2/PHD1), oxygen sensors mediating HIF hydroxylation. Validated for WB and ELISA; synthetic peptide immunogen, shipped cold. Ideal for hypoxia signaling studies.

[Learn More](#)

Parameter	Value
Product Name	EGLN1/EGLN2 Rabbit pAb
Remarks/Alias	EGLN1/EGLN2
Species	Human
GeneID (Human)	54583/112398
GeneID	54583/112398
Immunogen	A synthetic peptide corresponding to amino acids 250-350 of human EGLN1/EGLN2 (NP_071334.1/NP_444274.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9GZT9/Q96KS0
Protein Weight	36kDa/43kDa/46kDa/40kDa
Shipping	Ice bag