


Anti-Egr1 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P18146

Item Number: CM0006859



Introduction

Rabbit monoclonal antibody targeting human EGR1, validated for WB, IP, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for investigating transcriptional regulation, cell proliferation, and stress responses.

[Learn More](#)

Parameter	Value
Product Name	EGR1 Rabbit mAb
Remarks/Alias	TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268; EGR1
Species	Human
GeneID (Human)	1958
GeneID	1958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-250aa of human EGR1 (NP_001955.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18146
Protein Weight	58kDa
Shipping	Ice bag