


Anti-Egr1 Rabbit Polyclonal Antibody For Wb And Elisa - P18146

Item Number: CM0030291



Introduction

Rabbit polyclonal antibody against human EGR1, validated for WB and ELISA with cross-reactivity to mouse and rat. Suitable for studying transcriptional regulation, tumor suppression, and ischemia response.

[Learn More](#)

Parameter	Value
Product Name	EGR1 Rabbit pAb
Remarks/Alias	TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268; EGR1
Species	Human
GeneID (Human)	1958
GeneID	1958
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18146
Protein Weight	58 kDa
Shipping	Ice bag