


Anti-Egr1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P18146

Item Number: CM0028285



Introduction

Anti-EGR1 rabbit polyclonal antibody for ELISA, WB, IF/ICC. Detects human, mouse, rat EGR1 protein at 58kDa. Research tool for transcription regulation, cell growth, and stress response studies.

[Learn More](#)

Parameter	Value
Product Name	EGR1 Rabbit pAb
Remarks/Alias	TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268; EGR1
Species	Human
GeneID	1958
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EGR1 (NP_001955.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P18146
Protein Weight	58kDa
Shipping	Ice bag