

Anti-Egr2 Monoclonal Antibody For Wb, Ip, Elisa - P11161

Item Number: CM0005279

Introduction

EGR2 Rabbit monoclonal antibody validated for WB, IP, ELISA with cross-reactivity to human, mouse, and rat. Ideal for transcription and neurobiology research. - P11161

[Learn More](#)

Parameter	Value
Product Name	EGR2 Rabbit mAb
Remarks/Alias	AT591; CMT1D; CMT4E; KROX20; EGR2
Species	Human
Gene ID (Human)	1959
Gene ID	1959
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EGR2 (NP_000390.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11161
Protein Weight	50kDa
Shipping	Ice bag