

Anti-Egr3 Rabbit Polyclonal Antibody For Wb And Elisa - Q06889

Item Number: CM0029895

Introduction

Rabbit polyclonal antibody against human EGR3, validated for Western blot and ELISA. Cross-reacts with mouse. EGR3 is a zinc finger transcription factor involved in muscle spindle development. Suitable for research on neuromuscular diseases and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	EGR3 Rabbit pAb
Remarks/Alias	EGR-3; PILOT; EGR3
Species	Human
GeneID (Human)	1960
GeneID	1960
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EGR3 (NP_004421.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q06889
Protein Weight	43kDa
Shipping	Ice bag