

Anti-Egr4 Rabbit Polyclonal Antibody For Wb, Elisa - Q05215

Item Number: CM0029487

Introduction

Rabbit polyclonal antibody targeting human EGR4 (Early growth response protein 4) for Western Blot and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	EGR4 Rabbit pAb
Remarks/Alias	AT133; NGFIC; NGFI-C; PAT133; EGR4
Species	Human
GeneID (Human)	1961
GeneID	1961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-589 of human EGR4 (NP_001956.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05215
Protein Weight	62kDa
Shipping	Ice bag