

Anti-Egr4 Rabbit Polyclonal Antibody For Wb, Elisa - Q05215

Item Number: CM0018933

Introduction

Rabbit polyclonal antibody targeting Egr4 (also known as NGFIC, pAT133), validated for Western blotting and ELISA with mouse samples. Ideal for research on Egr4-mediated mitogenesis and differentiation.

[Learn More](#)

Parameter	Value
Product Name	Egr4 Rabbit pAb
Remarks/Alias	NGFIC; NGF1-C; NGFI-C; pAT133; Egr4
Species	Mouse
Gene ID (Human)	1961
Gene ID	13656
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-367 of mouse Egr4 (NP_065621.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q05215
Protein Weight	62kDa
Shipping	Ice bag