


Anti-Ehd1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H4M9

Item Number: CM0008639



Introduction

Anti-EHD1 Rabbit mAb validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat EHD1 (61kDa). Ideal for membrane trafficking, cilium biogenesis, and endocytic pathway research.

[Learn More](#)

Parameter	Value
Product Name	EHD1 Rabbit mAb
Remarks/Alias	PAST; PAST1; H-PAST; HPAST1; EHD1
Species	Human
GeneID (Human)	10938
GeneID	10938
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human EHD1 (Q9H4M9).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4M9
Protein Weight	61kDa
Shipping	Ice bag