

Anti-Ehd2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nzn4

Item Number: CM0015419

Introduction

Rabbit polyclonal antibody against human EHD2 (PAST2), an EH domain-containing protein. Validated for WB and ELISA in human, mouse, and rat samples. Supports studies on membrane trafficking and caveolae regulation.

[Learn More](#)

Parameter	Value
Product Name	EHD2 Rabbit pAb
Remarks/Alias	PAST2; EHD2
Species	Human
Gene ID (Human)	30846
Gene ID	30846
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 244-543 of human EHD2 (NP_055416.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NZN4
Protein Weight	61kDa
Shipping	Ice bag