

Anti-Ehf Polyclonal Antibody For Wb, Elisa - Q9Nzc4

Item Number: CM0027267

Introduction

Rabbit polyclonal antibody against human EHF, validated for Western blotting (WB) and ELISA. Detects human and mouse EHF; recombinant immunogen comprising amino acids 1-113. Suitable for research on epithelial differentiation and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	EHF Rabbit pAb
Remarks/Alias	ESE3; ESEJ; ESE3B; EHF
Species	Human
Gene ID	26298 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human EHF (NP_036285.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NZC4
Protein Weight	35kDa
Shipping	Ice bag