

Anti-Ehf Rabbit Polyclonal Antibody For Wb And Elisa - Q9Nzc4

Item Number: CM0016615

Introduction

Rabbit polyclonal antibody against human ETS homologous factor (EHF) validated for Western blot and ELISA. Cross-reacts with mouse. Useful for studying epithelial differentiation, transcriptional regulation, and tumor suppression mechanisms.

[Learn More](#)

Parameter	Value
Product Name	EHF Rabbit pAb
Remarks/Alias	ESE3, ESEJ, ESE3B, EHF
Species	Human
Gene ID (Human)	26298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human EHF (NP_036285.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NZC4
Protein Weight	35kDa
Shipping	Ice bag