

Anti-Ehhadh Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q08426

Item Number: CM0010062

Introduction

Rabbit monoclonal antibody against human EHHADH (peroxisomal bifunctional enzyme). Validated for WB, IF-P, IHC-P, ELISA; cross-reacts with mouse and rat. Ideal for fatty acid metabolism research.

[Learn More](#)

Parameter	Value
Product Name	EHHADH Rabbit mAb
Remarks/Alias	LBP; ECHD; LBFP; MFE1; PBFE; FRTS3; L-PBE
Species	Human
GeneID (Human)	1962
GeneID	1962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-284 of human EHHADH (NP_001957.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08426
Protein Weight	79kDa
Shipping	Ice bag