

Anti-Ehmt1 Polyclonal Antibody For Wb And Elisa - Q9H9B1

Item Number: CM0028280



Introduction

Rabbit polyclonal antibody against EHMT1 (GLP), validated for Western blot and ELISA. Detects human, mouse, and rat EHMT1 (~141 kDa). Suitable for chromatin biology and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	EHMT1 Rabbit pAb
Remarks/Alias	GLP; GLP1; KMT1D; KLEFS1; FP13812; EHMT1-IT1; EUHMTASE1; Eu-HMTase1; EHMT1
Species	Human
Gene ID (Human)	79813
Gene ID	79813
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human EHMT1 (NP_079033.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9B1
Protein Weight	141kDa
Shipping	Ice bag