

Anti-Ehmt2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96Kq7

Item Number: CM0003914



Introduction

Rabbit monoclonal antibody targeting EHMT2 (G9a) histone methyltransferase. Validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat EHMT2. Suitable for epigenetic transcriptional repression studies. - Q96KQ7

[Learn More](#)

Parameter	Value
Product Name	EHMT2 Rabbit mAb
Remarks/Alias	G9A; BAT8; GAT8; NG36; KMT1C; C6orf30; EHMT2
Species	Human
GeneID (Human)	10919
GeneID	10919
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-300 of human EHMT2 (Q96KQ7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96KQ7
Protein Weight	132 kDa
Shipping	Ice bag