

Anti-Ei24 Rabbit Polyclonal Antibody For Wb, Elisa - O14681

Item Number: CM0022056

Introduction

Rabbit polyclonal antibody targeting human EI24 (EPG4/PIG8), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for studying p53-mediated apoptosis and autophagy pathways.

[Learn More](#)

Parameter	Value
Product Name	EI24 Rabbit pAb
Remarks/Alias	EPG4; PIG8; TP53I8; EI24
Species	Human
Gene ID (Human)	9538
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 240-340 of human EI24 (NP_004870.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14681
Protein Weight	39kDa
Shipping	Ice bag