


Anti-Eif1A Monoclonal Antibody For Wb, Elisa - P47813

Item Number: CM0003009



Introduction

Rabbit monoclonal antibody targeting human eIF1A (EIF1AX) for WB and ELISA applications. Cross-reacts with mouse. Key translation initiation factor involved in 43S complex scanning.

[Learn More](#)

Parameter	Value
Product Name	eIF1A Rabbit mAb
Remarks/Alias	EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1; eIF1A
Species	Human
GeneID (Human)	1964
GeneID	1964
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 65-144 of human eIF1A (P47813).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47813
Protein Weight	16kDa
Shipping	Ice bag