


Anti-Eif1Ax Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P47813

Item Number: CM0027169



Introduction

Rabbit polyclonal EIF1AX antibody (P47813) validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat eIF1AX, a 43S pre-initiation complex component essential for translation initiation. Ideal for mechanistic studies.

[Learn More](#)

Parameter	Details
Product Name	EIF1AX Rabbit pAb
Remarks/Alias	EIF1A, EIF4C, eIF-1A, eIF-4C, EIF1AP1, EIF1AX
Species	Human
GeneID (Human)	1964
GeneID	1964
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-144 of human EIF1AX (NP_001403.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47813
Protein Weight	16kDa
Shipping	Ice bag