

## Anti-Eif1Ay Rabbit Polyclonal Antibody For Wb, Elisa - O14602

Item Number: CM0025448

### Introduction

Rabbit polyclonal antibody targeting human EIF1AY (eIF-4C). Validated for WB and ELISA, with cross-reactivity to mouse and rat, suitable for translation initiation research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | EIF1AY Rabbit pAb                                                                                                  |
| Remarks/Alias    | eIF-4C; EIF1AY                                                                                                     |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 9086                                                                                                               |
| GeneID           | 9086                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-144 of human EIF1AY (NP_004672.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | O14602                                                                                                             |
| Protein Weight   | 16kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |