


Anti-Eif2A Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P05198

Item Number: CM0007401



Introduction

High-quality Anti-eIF2 α Rabbit Monoclonal Antibody validated for WB, IF/ICC, IHC-P, ELISA applications in Human, Mouse, Rat samples. Explore eukaryotic translation initiation factor 2 subunit 1 function in stress response and protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	eIF2 α Rabbit mAb
Remarks/Alias	EIF2; EIF-2; EIF2A; EIF-2A; EIF-2alpha; eIF2 α
Species	Human
GeneID (Human)	1965
GeneID	1965
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-315 of human eIF2 α (NP_004085.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05198
Protein Weight	36kDa
Shipping	Ice bag