


Anti-Eif2A Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9By44

Item Number: CM0010606



Introduction

Rabbit monoclonal antibody against eukaryotic translation initiation factor 2A (eIF2A), validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Suitable for translation research. MW 65 kDa.

[Learn More](#)

Parameter	Value
Product Name	eIF2A Rabbit mAb
Remarks/Alias	CDA02; EIF-2A; MST089; MSTP004; MSTP089; eIF2A
Species	Human
GeneID	83939
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 450-550 of human eIF2A (Q9BY44)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BY44
Protein Weight	65kDa
Shipping	Ice bag