


Anti-Eif2A Rabbit Polyclonal Antibody For Wb, Elisa - Q9By44

Item Number: CM0020240



Introduction

CamelBio's rabbit polyclonal anti-eIF2A antibody is validated for WB and ELISA, cross-reacting with human and mouse. Supports translation initiation and protein synthesis research.

[Learn More](#)

Parameter	Value
Product Name	eIF2A Rabbit pAb
Remarks/Alias	CDA02; EIF-2A; MST089; MSTP004; MSTP089; eIF2A
Species	Human
GeneID (human)	83939
GeneID	83939
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF2A (NP_114414.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BY44
Protein Weight	65kDa
Shipping	Ice bag