


## Anti-Eif2A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P05198

Item Number: CM0015227



### Introduction

Rabbit polyclonal antibody against human eIF2 $\alpha$  (eukaryotic translation initiation factor 2 subunit alpha). Validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying translation initiation, integrated stress response, and mitophagy.

[Learn More](#)

| Parameter        | Value                                                                                                                                          |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | eIF2 $\alpha$ Rabbit pAb                                                                                                                       |
| Remarks/Alias    | EIF2, EIF-2, EIF2A, EIF-2A, EIF-2alpha, eIF2 $\alpha$                                                                                          |
| Species          | Human                                                                                                                                          |
| GeneID (Human)   | 1965                                                                                                                                           |
| GeneID           | 1965                                                                                                                                           |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human eIF2 $\alpha$ (NP_004085.1). |
| Source           | Rabbit                                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                                          |
| Application      | WB, IHC-P, ELISA                                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                              |
| SWISS            | P05198                                                                                                                                         |
| Protein Weight   | 36 kDa                                                                                                                                         |
| Shipping         | Ice bag                                                                                                                                        |