

Anti-Eif2Ak1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bqi3

Item Number: CM0016582

Introduction

Rabbit polyclonal antibody targeting human EIF2AK1 (HRI), validated for WB, IHC-P, ELISA; cross-reacts with rat. Useful for integrated stress response, heme deficiency, and mitophagy research.

[Learn More](#)

Parameter	Value
Product Name	EIF2AK1 Rabbit pAb
Remarks/Alias	HCR; HRI; hHRI; LEMSPAD; EIF2AK1
Species	Human
Gene ID	27102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 337-586 of human EIF2AK1 (NP_055228.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BQI3
Protein Weight	71kDa
Shipping	Ice bag