

Anti-Eif2Ak1 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Bqi3

Item Number: CM0007031

Introduction

Rabbit monoclonal antibody against human EIF2AK1, validated for Western blot and ELISA. Targets heme-regulated inhibitor kinase involved in stress responses, hemoglobin synthesis, and mitochondrial dysfunction. Suitable for research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] EIF2AK1 Rabbit mAb
Remarks/Alias	HCR; HRI; hHRI; LEMSPAD; [KD Validated] EIF2AK1
Species	Human
Gene ID (Human)	27102
Gene ID	27102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 337-586 of human EIF2AK1 (NP_055228.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BQI3
Protein Weight	71kDa
Shipping	Ice bag