

Anti-Eif2Ak1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bqi3

Item Number: CM0027860

Introduction

High-quality rabbit polyclonal antibody against EIF2AK1 (HRI) for WB and ELISA applications. Reactive with human and mouse samples. Ideal for studying integrated stress response, heme sensing, and translational control.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] EIF2AK1 Rabbit pAb
Remarks/Alias	HCR; HRI; hHRI; LEMSPAD; K1
Species	Human
Gene ID (Human)	27102
Gene ID	27102
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 371-630 of human EIF2AK1 (NP_055228.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BQI3
Protein Weight	71kDa
Shipping	Ice bag