


## Anti-Eif2B Epsilon (Eif2B5) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13144

Item Number: CM0010409



### Introduction

Rabbit monoclonal antibody targeting human eIF2B epsilon (EIF2B5), validated for WB, IHC-P, IF/ICC, ELISA. Essential for studying translation initiation and stress responses. UniProt Q13144.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | eIF2B epsilon (EIF2B5) Rabbit mAb                                                                               |
| Remarks/Alias    | CLE; CACH; LVWM; EIF-2B; EIF2Bepsilon; eIF2B epsilon (EIF2B5)                                                   |
| Species          | Human                                                                                                           |
| Gene ID (Human)  | 8893                                                                                                            |
| Gene ID          | 8893                                                                                                            |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-100 of human eIF2B epsilon (EIF2B5) (Q13144) |
| Source           | Rabbit                                                                                                          |
| Category         | Monoclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                        |
| Cross Reactivity | Human                                                                                                           |
| SWISS            | Q13144                                                                                                          |
| Protein Weight   | 80kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |