


Anti-Eif2B Epsilon (Eif2B5) Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q13144

Item Number: CM0017839



Introduction

Rabbit polyclonal antibody against human eIF2B epsilon (EIF2B5). Suitable for WB, IF/ICC, IP, ELISA; cross-reacts with human, mouse, rat. Essential for studying translation initiation and eIF2B complex function.

[Learn More](#)

Product Name	eIF2B epsilon (EIF2B5) Rabbit pAb
Remarks/Alias	CLE; CACH; LVWM; EIF-2B; EIF2Bepsilon; eIF2B epsilon (EIF2B5)
Species	Human
Gene ID (Human)	8893
Gene ID	8893
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 442-721 of human eIF2B epsilon (EIF2B5) (NP_003898.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13144
Protein Weight	80kDa
Shipping	Ice bag