


Anti-Eif2B1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q14232

Item Number: CM0029132



Introduction

Rabbit polyclonal antibody targeting human EIF2B1, validated for Western blot, IF/ICC, and ELISA, with cross-reactivity to mouse. Ideal for studying translation regulation, eIF2B complex function, and leukodystrophy research.

[Learn More](#)

Parameter	Value
Product Name	EIF2B1 Rabbit pAb
Remarks/Alias	EIF2B; EIF2BA; EIF2Balpha; EIF2B1
Species	Human
Gene ID (Human)	1967
Gene ID	1967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-305 of human EIF2B1 (NP_001405.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14232
Protein Weight	34kDa
Shipping	Ice bag