


Anti-Eif2B2/Eif2Bbeta Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49770

Item Number: CM0029937



Introduction

Rabbit polyclonal antibody targeting human EIF2B2/eIF2Bbeta, validated for WB, IHC-P, IF/ICC, IP, and ELISA applications. Cross-reacts with mouse and rat. Ideal for translation initiation research.

[Learn More](#)

Parameter	Value
Product Name	EIF2B2 Rabbit pAb
Remarks/Alias	EIF2B; EIF2Bbeta; EIF-2Bbeta; EIF2B2
Species	Human
GeneID (Human)	8892
GeneID	8892
Immunogen	Recombinant fusion protein containing amino acids 1-120 of human EIF2B2 (NP_055054.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49770
Protein Weight	39kDa
Shipping	Ice bag