


Anti-Eif2B3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nr50

Item Number: CM0017255



Introduction

Rabbit polyclonal antibody targeting EIF2B3, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Recognizes the eIF2B subunit gamma, crucial for translation initiation and stress response. Ideal for translation studies.

[Learn More](#)

Parameter	Value
Product Name	EIF2B3 Rabbit pAb
Remarks/Alias	EIF-2B; EIF2Bgamma; EIF2B3
Species	Human
GeneID (human)	8891
GeneID	8891
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human EIF2B3 (NP_065098.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NR50
Protein Weight	50kDa
Shipping	Ice bag