


Anti-Eif2B4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9UI10

Item Number: CM0022835



Introduction

Anti-EIF2B4 Rabbit Polyclonal Antibody (Q9UI10) validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Ideal for investigating translation initiation factor eIF2B complex function and regulation.

[Learn More](#)

Parameter	Value
Product Name	EIF2B4 Rabbit pAb
Remarks/Alias	EIF2B; EIF-2B; EIF2Bdelta; EIF2B4
Species	Human
GeneID (Human)	8890
GeneID	8890
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-100 of human EIF2B4 (Q9UI10)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UI10
Protein Weight	58kDa
Shipping	Ice bag