

Anti-Eif2C3 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9H9G7

Item Number: CM0003690



Introduction

Rabbit monoclonal antibody against human eIF2C3 (Argonaute-3), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Designed for RNAi and gene silencing research applications.

[Learn More](#)

Parameter	Value
Product Name	eIF2C3 Rabbit mAb
Remarks/Alias	EIF2C3; Argonaute3; hAgo3; eIF2C3
Species	Human
Gene ID (Human)	192669
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 721-860 of human eIF2C3 (Q9H9G7)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H9G7
Protein Weight	97kDa
Shipping	Ice bag