

Anti-Eif2S3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P41091

Item Number: CM0025198



Introduction

Rabbit polyclonal antibody targeting human EIF2S3 (eIF2 γ), validated for Western blot, IHC-P, and ELISA. Demonstrates cross-reactivity with mouse and rat. Ideal for investigating translation regulation and eIF2 function.

[Learn More](#)

Parameter	Value
Product Name	EIF2S3 Rabbit pAb
Remarks/Alias	EIF2; EIF2G; MEHMO; MRXSBRK; eIF-2gA; EIF2gamma; EIF2S3
Species	Human
Gene ID (Human)	1968
Gene ID	1968
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human EIF2S3 (NP_001406.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41091
Protein Weight	51kDa
Shipping	Ice bag